

EXPERIMENT PROCEDURE

- Determine the two positions of a thin lens where a sharp image is formed.
- Determine the focal length of a thin lens.

OBJECTIVE

Determine the focal length of a lens using the Bessel method

SUMMARY

On an optical bench it is possible to set up a light source, a lens, a screen and an object to be imaged in such a way that a well focussed image appears on the screen. Using the geometric relationships between the ray paths for a thin lens, it is possible to determine its focal length.

REQUIRED APPARATUS

| Quantity | Description | Number |
|----------|---|------------|
| 1 | Optical Bench K, 1000 mm | 1009696 |
| 4 | Optical Rider K | 1000862 |
| 1 | Optical Lamp K | 1000863 |
| 1 | Transformer 12 V, 25 VA (230 V, 50/60 Hz) | 1000866 or |
| | Transformer 12 V, 25 VA (115 V, 50/60 Hz) | 1000865 |
| 1 | Convex Lens K, $f = 50$ mm | 1000869 |
| 1 | Convex Lens K, $f = 100$ mm | 1010300 |
| 1 | Clamp K | 1008518 |
| 1 | Set of 4 Image Objects | 1000886 |
| 1 | Projection Screen K, White | 1000879 |

1

BASIC PRINCIPLES

The focal length f of a lens refers to the distance between the main plane of the lens and its focal point, see Fig.1. This can be determined using the Bessel method (devised by Friedrich Wilhelm Bessel). This involves measuring the various separations between the optical components on the optical bench.

From Fig.1 and Fig. 2 it can be seen that the following relationship must apply for a thin lens:

$$(1) \quad a = b + g$$

a : distance between object G and image B
 b : distance between lens and image B
 g : distance between object G and lens

By plugging these values into the lens equation

$$(2) \quad \frac{1}{f} = \frac{1}{b} + \frac{1}{g}$$

f : focal length of lens

the following is obtained:

$$(3) \quad \frac{1}{f} = \frac{a}{a \cdot g - g^2}$$

This corresponds to a quadratic equation with the following pair of solutions:

$$(4) \quad g_{1,2} = \frac{a}{2} \pm \sqrt{\frac{a^2}{4} - a \cdot f}$$

A sharp image is obtained for each of the object distances g_1 and g_2 . The difference e between them allows the focal length to be determined:

$$(5) \quad e = g_1 - g_2 = \sqrt{a^2 - 4af}$$

The difference e is the difference between the two lens positions P_1 and P_2 , which result in a focussed image.

EVALUATION

A formula for the focal length of a thin lens can be derived using the

Bessel method from equation (4) $f = \frac{a^2 - e^2}{4a}$




Fig. 1: Schematic showing the definition of focal length for a thin lens




Fig.2: Schematic of ray paths through a lens




Fig.3: Schematic showing the two lens positions which result in a well focussed image on the screen